

Anti-ADAMTS2 antibody



Catalog Number: 176458

Product name

Anti-ADAMTS2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ADAMTS2

Preparation

This antigen of this antibody was synthesized peptide derived from adamts-2 . at aa range: 1140-1220

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of HepG2 K562 HeLa Colo using ADAMTS-2 antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000