

Anti-Monocarboxylic acid transporter 1 Rabbit antibody



Catalog Number: 168858

Product name

Monoclonal

Anti-Monocarboxylic acid transporter 1 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

Monocarboxylic acid transporter 1 Rabbit monoclonal antibody


Dilutions

1:1000

Preparation

Antigen: A synthetic peptide of human Monocarboxylic acid transporter 1

Validations



Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Monocarboxylic acid transporter 1 in K562,HeLa cell lysates using Monocarboxylic acid transporter 1 Rabbit mAb(1:1000 diluted).Predicted band size:54kDa.Observed band size:54kDa.

Clonality