

Human IDI2 (His tag) recombinant protein



Catalog Number: 518858

General Information

Gene Name Synonym

Isopentenyl-Diphosphate Delta-Isomerase 2, Isopentenyl Pyrophosphate Isomerase 2, IPP Isomerase 2, IPPI2, IDI2

Protein Construction

A DNA sequence encoding the Human IDI2 Met1-Val227 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, 0.1mM PMSF, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human IDI2 (His tag) recombinant protein