

# SOCS1 Polyclonal Antibody



Catalog Number: 160158

## Product name

SOCS1 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to SOCS1

## Preparation

Antigen: A synthetic peptide of human SOCS1

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF, FC

## Dilutions


WB 1:500 - 1:1000

IHC 1:50 - 1:200

IF 1:20 - 1:50

FC 1:20 - 1:50

## Validations



## Western blot - SOCS1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SOCS1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.