

Anti-RORB antibody



Catalog Number: 142158

Product name

Anti-RORB antibody

Specificity

Human RORB

Antibody description

Rabbit Polyclonal to Human RORB

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human RORB, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 10-30 µg/ml

Validations



Items	Lanes	A
Sample (whole cell lysate)		293T
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10-30 µg/ml
Secondary Antibody		Dylight 800-labeled Antibody to Rabbit IgG (1:14), at 1:5000 dilution. Developed using Odyssey imaging system.
Explanation		Predicted band size : 53 kDa Observed band size : 58 kDa

RORB Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot