

# BMP1 Polyclonal Antibody



Catalog Number: 163158

---

## Product name

BMP1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to BMP1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 725-986 of human BMP1 (NP\_006120.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - BMP1 Polyclonal Antibody  
Immunofluorescence analysis of HeLa cells using BMP1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.