

p63 antibody



Catalog Number: 113558

Product name

p63 antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:200-1:1000

IHC: 1:20-1:200

Antibody description

p63 Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive IHC detected in human tonsillitis tissue, human breast hyperplasia tissue, human lung cancer tissue. Positive WB detected in A431 cells, Apoptosis HeLa cells, human heart tissue, mouse heart tissue, U-937 cells. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 63-68 kDa

IF: 1:50-1:500

Validations



Preparation

This antibody was obtained by immunization of p63 recombinant protein (Accession Number: NM_003722). Purification method: Antigen affinity purified.

A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113558 (TP63 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF


Dilutions




IP Result of anti-TP63 (IP: Catalog No: 113558, 4ug; Detection: Catalog No: 113558 1:300) with A431 cells lysate 1200ug.

p63 antibody


Catalog Number: 113558



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:113558(TP63 Antibody) at dilution of 1:200 (under 40x lens)



Immunofluorescent analysis of (10% Formaldehyde) fixed A431 cells using Catalog No:113558(TP63 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human breast hyperplasia tissue slide using Catalog No:113558(TP63 Antibody) at dilution of 1:200 (under 40x lens)