

Anti-Adiponectin Rabbit antibody



Catalog Number: 168558

Product name

Anti-Adiponectin Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IP

Antibody description

Adiponectin Rabbit monoclonal antibody

Dilutions

1:1000-1:5000

Preparation


Antigen: A synthetic peptide of human Adiponectin

1:20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Adiponectin in K562, HeLa cell lysates using Adiponectin Rabbit mAb(1:1000 diluted). Predicted band size: 26kDa. Observed band size: 30kDa.

Clonality

Monoclonal