

Human THOP1 (His tag) recombinant protein



Catalog Number: 519558

General Information

Gene Name Synonym

Thimet Oligopeptidase, Endopeptidase 24.15, MP78, THOP1

Protein Construction

A DNA sequence encoding the Human THOP1 Lys2-Cys689 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 100mM NaCl, 10% Glycerol, pH 7.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human THOP1 (His tag) recombinant protein