

# GGCT antibody

Catalog Number: 110958

## Product name

GGCT antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GGCT Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse kidney tissue, mouse liver tissue, U-937 cells. Positive IHC detected in human liver cancer. Observed molecular weight by Western-blot: 21-25kd

## Preparation

This antibody was obtained by immunization of GGCT recombinant protein (Accession Number: NM\_024051). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110958(GGCT antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:110958(GGCT antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:110958(GGCT antibody) at dilution of 1:50

# GGCT antibody



Catalog Number: 110958

---

(under 40x lens)