

# PLEKHA8 antibody



Catalog Number: 113958

## Product name

PLEKHA8 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PLEKHA8 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IHC detected in human kidney. Observed molecular weight by Western-blot: 49kd

## Preparation

This antibody was obtained by immunization of PLEKHA8 recombinant protein (Accession Number: NM\_032639). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113958(PLEKHA8 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113958(PLEKHA8 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113958(PLEKHA8 antibody) at dilution of 1:50 (under 40x lens)