

ZNF2 antibody

Catalog Number: 116958

Product name

ZNF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZNF2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in human spleen tissue, mouse spleen tissue. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of ZNF2 recombinant protein (Accession Number: NM_001017396). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116958(ZNF2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116958(ZNF2 antibody) at dilution of 1:50 (under 10x lens)