

Anti-PIP5K1C antibody



Catalog Number: 176258

Product name

Anti-PIP5K1C antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PIP5K1C

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human pip5k1c. aa range:305-354

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using PIPK I γ Polyclonal Antibody