

# Anti-DR5 (7F4) Mouse antibody



Catalog Number: 168258

## Product name

Anti-DR5 (7F4) Mouse antibody

## Specificity

Human, Mouse

## Antibody description

DR5 (7F4) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human DR5 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF

## Dilutions


WB: 1/500-1/2000

;ICC: 1/100

## Validations



Western blot detection of DR5 in Mouse kidney, Mouse brain, HCT116 and HeLa cell lysates using DR5 mouse mAb (1:500-1:2000 diluted). Predicted band size: 40/48kDa. Observed



band size: 40/48kDa.

Immunocytochemistry of HeLa cells fixed by Paraformaldehyde and using DR5 mouse mAb diluted 1:100.