

CENPP antibody



Catalog Number: 109258

Product name

CENPP antibody

Immunogen

[Human CENPP Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

CENPP Rabbit Polyclonal antibody. Positive WB detected in HEK-293T cells. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of CENPP recombinant protein (Accession Number: NM_001012267). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HEK-293T cells were subjected to SDS PAGE followed by western blot with Catalog No:109258(CENPP Antibody) at dilution of 1:600