

# LTB4R antibody

Catalog Number: 112358

---

## Product name

LTB4R antibody

## Specificity

Human; other species not tested.

## Antibody description

LTB4R Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human breast cancer tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_181657).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112358(LTB4R Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112358(LTB4R Antibody) at dilution of 1:50 (under 40x lens)