

# NR5A2 Polyclonal Antibody



Catalog Number: 163358

---

## Product name

NR5A2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to NR5A2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human NR5A2 (NP\_003813.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions

IHC 1:50 - 1:200

## Validations



Immunohistochemistry - NR5A2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using NR5A2 antibody at dilution of 1:100 (40x lens).