

# Anti-REG1B antibody



Catalog Number: 106358

## Product name

Anti-REG1B antibody

## Immunogen

[Human REG1B \(His Tag\) recombinant protein](#)

## Specificity

Human REG1B

## Antibody description

Rabbit monoclonal to REG1B

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human REG1B (rh REG1B; NP\_006498.1; Met1-Asn166).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-10 µg/mL

## Validations



REG1B Antibody, Rabbit MAb,  
Immunohistochemistry

Immunochemical staining of human REG1B in human pancreas with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).