

TRERF1 antibody



Catalog Number: 116358

Product name

TRERF1 antibody

Specificity

human, mouse; other species not tested.

Antibody description

TRERF1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 132-140kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_033502). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116358(TRERF1 antibody) at dilution of 1:300