

# Human RAET1E / ULBP4 (Fc Tag) recombinant protein



Catalog Number: 501758

## General Information

### Gene Name Synonym

Lymphocyte effector toxicity activation ligand; RAE-1-like transcript 4; Retinoic acid early transcript 1E

### Protein Construction

A DNA sequence encoding the human RAET1E (NP\_631904.1) (Met1-Asp225) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

His 31

### Molecular Mass

The recombinant human RAET1E consists 433 amino acids and predicts a molecular mass of 49.1 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human RAET1E / ULBP4 Protein (Fc Tag) SDS-PAGE