

Visfatin antibody

Catalog Number: 116758

Product name

Visfatin antibody

Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

Antibody description

Visfatin Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue, HL-60 cells, mouse liver tissue, rat liver tissue, zebrafish tissue. Observed molecular weight by Western-blot: 46-52 kDa

Preparation

This antibody was obtained by immunization of Visfatin recombinant protein (Accession Number: BC020691). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:20-1:200

Validations




Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116758(NAMPT antibody) at dilution of 1:50



(under 10x lens)

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116758(NAMPT antibody) at dilution of 1:1000



IP Result of anti-NAMPT (IP:Catalog No:116758, 3ug; Detection:Catalog No:116758 1:1000) with mouse heart tissue lysate 4000ug.