

# Anti-HJURP Rabbit antibody



Catalog Number: 167758

## Product name

Anti-HJURP Rabbit antibody

## Specificity

Human

## Antibody description

HJURP Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human HJURP

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/1000

IHC: 1/20

## Validations



Western blot detection of HJURP in HeLa,CHO-K1,  
using HJURP Rabbit mAb(1:1000 diluted)