

# DPP3 antibody



Catalog Number: 110068

## Product name

DPP3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DPP3 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, human placenta tissue. Observed molecular weight by Western-blot: 69 kDa, 83 kDa

## Preparation

This antibody was obtained by immunization of DPP3 recombinant protein (Accession Number: NM\_130443). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:110068(DPP3 antibody) at dilution of 1:1500