

# ApoE Polyclonal Antibody



Catalog Number: 161068

---

## Product name

ApoE Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ApoE

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ApoE (NP\_000032.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions

IHC 1:50 - 1:100

## Validations



Immunohistochemistry - ApoE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using ApoE antibody at dilution of 1:100 (40x lens).