

TLE1 antibody

Catalog Number: 116068

Product name

TLE1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TLE1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue.

Preparation

This antibody was obtained by immunization of TLE1 recombinant protein (Accession Number: NM_005077). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116068(TLE1 antibody) at dilution of 1:50 (under 10x lens)