

Anti-BTK (5B12) Mouse antibody



Catalog Number: 168068

Product name

Anti-BTK (5B12) Mouse antibody

Specificity

Human

Antibody description

BTK (5B12) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human BTK protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB, IP


Dilutions

WB: 1/1000

Validations



Western blot detection of BTK in Daudi, Jurkat (BTK negative), K562, Ramos and Raji cell lysates using BTK mouse mAb (1:1000 diluted). Predicted band size: 77KDa. Observed



band size: 77KDa.

Immunoprecipitation analysis of Raji cell lysates using BTK mouse mAb.