

Anti-BTK (5B12) Mouse antibody



Catalog Number: 168068

Product name

WB, IP

Anti-BTK (5B12) Mouse antibody

Dilutions

Specificity

WB: 1/1000

Human

Antibody description

BTK (5B12) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human BTK protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Western blot detection of BTK in Daudi, Jurkat (BTK negative), K562, Ramos and Raji cell lysates using BTK mouse mAb (1:1000 diluted). Predicted band size: 77KDa. Observed



band size: 77KDa.

Immunoprecipitation analysis of Raji cell lysates using BTK mouse mAb.