

Anti-PNLIP antibody

Catalog Number: 100468

Product name

Anti-PNLIP antibody

Immunogen

[Human PNLIP \(His Tag\) recombinant protein](#)

Specificity

Human PNLIP

Antibody description

Rabbit monoclonal to PNLIP

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human PNLIP (rh PNLIP; P16233; Met1-Cys465).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.5-5 µg/mL

Validations



PNLIP Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of human PNLIP in human pancreas with rabbit monoclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).