

Rat COPZ1/Zeta-1 COP (His Tag) recombinant protein



Catalog Number: 501468

General Information

Protein Construction

A DNA sequence encoding the rat COPZ1 (NP_057141.1) (Met1-Arg177) was expressed with a polyhistidine tag at the N-terminus.

Organism

Rat

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant rat COPZ1 consists 195 amino acids and predicts a molecular mass of 22.4 kDa.

Formulation

Lyophilized from sterile 50 mM Tris, 50 mM Arg, pH 8.5.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat COPZ1 Protein (His Tag)