

# TERF1 Polyclonal Antibody



Catalog Number: 164468

## Product name

TERF1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to TERF1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 240-419 of human TERF1 (NP\_003209.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Western blot - TERF1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TERF1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.