

Anti-IKK alpha (4C8) Mouse antibody



Catalog Number: 168468

Product name

Mouse IgG1

Anti-IKK alpha (4C8) Mouse antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description

WB: 1/500 - 1/2000

IKK alpha (4C8) Mouse monoclonal antibody

ELISA: 1/10000

Preparation

Antigen: Purified recombinant fragment of human CHUK expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot analysis of extracts from Jurkat and Raji cell lysates using IKK α mouse mAb (1:500 diluted). Predicted band size:85KDa. Observed band size:85KDa.