

Interferon alpha antibody



Catalog Number: 111868

Product name

Interferon alpha antibody

Specificity

Human; other species not tested.

Antibody description

Interferon alpha Rabbit Polyclonal antibody.
Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of Interferon alpha recombinant protein (Accession Number: NM_024013). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:111868(IFNA1 Antibody) at dilution of 1:3000