

# Rosopatamab tetraxetan (Synonyms: DOTA-HUJ-591-GS; DOTA-HUJ-591)



Catalog Number: 176868

## Product name

Rosopatamab tetraxetan (Synonyms: DOTA-HUJ-591-GS; DOTA-HUJ-591)

## Specificity

Human FOLH1 / PSMA

## Antibody description

Anti-FOLH1 / PSMA Antibody (Rosopatamab tetraxetan)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-FOLH1 / PSMA Antibody (rosopatamab tetraxetan) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-FOLH1 / PSMA Antibody (rosopatamab tetraxetan) is more than 95%, determined by SEC-HPLC.