

Anti-muscle Actin Rabbit antibody



Catalog Number: 167868

Product name

Anti-muscle Actin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

muscle Actin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human muscle Actin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions


WB: 1/1000

IHC: 1/201/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of muscle Actin in HeLa, CHO-K1, C6, Rat Brain lysates using muscle Actin Rabbit mAb(1:1000 diluted). Predicted band size:42kDa. Observed band size:42kDa.