

PRC1 antibody

Catalog Number: 114168

Product name

Polyclonal

PRC1 antibody

Ig Type**Specificity**

Rabbit IgG

Human,Mouse,Rat; other species not tested.

Applications**Antibody description**

ELISA, WB

PRC1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293T cells, HeLa cells, HepG2 cells, MCF-7 cells. Observed molecular weight by Western-blot: 67 kDa


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of PRC1 recombinant protein (Accession Number: BC005140). Purification method: Antigen affinity purified.

Validations**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HEK-293T cells were subjected to SDS PAGE followed by western blot with Catalog No:114168(PRC1 Antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality