

CDK9 antibody



Catalog Number: 109168

Product name

CDK9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CDK9 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in A431 cells, HEK-293 cells, human placenta tissue, Jurkat cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of CDK9 recombinant protein (Accession Number: NM_001261). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109168(CDK9 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109168(CDK9 antibody) at dilution of 1:50 (under 40x lens)