

# Anti-PLXNA1 antibody



Catalog Number: 176568

## Product name

Anti-PLXNA1 antibody

ICC: 1:50-1:200

FC: 1:50-1:100

## Specificity

Human, Mouse, Rat

WB: 1:500

## Antibody description

Rabbit monoclonal antibody to PLXNA1

## Preparation

This antigen of this antibody was recombinant protein

## Formulation

Liquid, 1\*TBS (pH7.4), 0.05% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC, FC

## Dilutions

## Validations




Fig1: ICC staining Plexin A1 in LOVO cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig2: Flow cytometric analysis of LOVO cells with Plexin A1 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody.