

Anti-Ataxin 3 Rabbit antibody



Catalog Number: 167568

Product name

Anti-Ataxin 3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Ataxin 3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Ataxin 3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Ataxin 3 in HeLa,A549,HL-60,U251,U87-MG cell lysates using Ataxin 3 Rabbit mAb(1:1000 diluted).Predicted band size:42kDa.Observed band size:42kDa.