

BPGM antibody



Catalog Number: 108568

Product name

BPGM antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BPGM Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human blood tissue, rat liver tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of BPGM recombinant protein (Accession Number: XM_011516527). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:117215(BPGM Antibody) at dilution of 1:1000