

GHRL antibody



Catalog Number: 110968

Product name

GHRL antibody

Immunogen

[Human GHRL Recombinant protein \(GST tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

GHRL Rabbit Polyclonal antibody. Positive WB detected in human stomach tissue. Positive IHC detected in human stomach tissue, mouse stomach tissue, rat stomach tissue. Observed molecular weight by Western-blot: 13 kDa

Preparation

This antibody was obtained by immunization of GHRL recombinant protein (Accession Number: NM_001302825). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:110968(GHRL antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human stomach using Catalog No:110968(GHRL antibody) at dilution of 1:100 (under 10x lens)



GHRL antibody



Catalog Number: 110968

Immunohistochemical of paraffin-embedded

human stomach using Catalog No:110968(GHRL antibody) at dilution of 1:100 (under 40x lens)