

KCNC1 Polyclonal Antibody



Catalog Number: 163968

Product name

KCNC1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to KCNC1

Preparation

Antigen: A synthetic peptide of human KCNC1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:25 - 1:100

Validations



Western blot - KCNC1 Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using KCNC1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.