

RAB3C-specific antibody



Catalog Number: 106968

Product name

RAB3C-specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB3C-specific Rabbit Polyclonal antibody.
Positive WB detected in human brain tissue.
Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_138453).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114443(RAB3C-specific antibody) at dilution of



1:300
Immunofluorescent analysis of HepG2 cells using Catalog No:114443(RAB3C-specific Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG