

NR3C2 antibody

Catalog Number: 113268

Product name

NR3C2 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

NR3C2 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, HeLa cells. Positive IHC detected in human kidney tissue, human colon tissue, human skin tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 94-100kd

Validations



Preparation

This antibody was obtained by immunization of NR3C2 recombinant protein (Accession Number: NM_001354819). Purification method: Antigen affinity purified.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113268(NR3C2 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:




Immunohistochemical of paraffin-embedded human kidney using Catalog No:113268(NR3C2 antibody) at dilution of 1:50 (under 10x lens)



NR3C2 antibody


Catalog Number: 113268

Immunohistochemical of paraffin-embedded human kidney using Catalog No:113268(NR3C2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using NR3C2 antibody Catalog No:113268 at 1:25

dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-NR3C2 (IP:Catalog No:113268, 4ug; Detection:Catalog No:113268 1:500) with HEK-293 cells lysate 3200ug.