

Anti-ANP32C antibody

Catalog Number: 176268

Product name

Anti-ANP32C antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ANP32C

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human anp32c. aa range:121-170

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using PP32R1 Polyclonal Antibody diluted at 1: 1000