

Anti-ULK3 Rabbit antibody



Catalog Number: 167268

Product name

Anti-ULK3 Rabbit antibody

Specificity

Human

Antibody description

ULK3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human ULK3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/20-1/50

FC: 1/20

Validations



Western blot detection of ULK3 in K562,HeLa cell lysates using ULK3 Rabbit pAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.