

# HSP90 antibody

Catalog Number: 107268

## Product name

HSP90 antibody

## Immunogen

[Human HSP90 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

HSP90 Mouse Monoclonal antibody. Positive WB detected in Y79 cells. Positive FC detected in K-562 cells. Observed molecular weight by Western-blot: 85-90 kDa

## Preparation

This antibody was obtained by immunization of HSP90 recombinant protein (Accession Number: BC023006). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

ELISA, WB, FC


## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:107268(HSP90 Antibody) at dilution of 1:4000



1X10<sup>6</sup> K-562 cells were stained with 0.2ug HSP90 antibody (Catalog No:107268, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.