

PDE7A antibody



Catalog Number: 113668

Product name

PDE7A antibody

Specificity

Human; other species not tested.

Antibody description

PDE7A Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of PDE7A recombinant protein (Accession Number: NM_002603). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113668(PDE7A Antibody) at dilution of 1:1000