

# A4GNT antibody



Catalog Number: 107668

## Product name

A4GNT antibody

## Immunogen

[Human A4GNT Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

A4GNT Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, L02 cells. Positive IHC detected in human stomach tissue, human pancreas tissue. Observed molecular weight by Western-blot: 37-39kd

## Preparation

This antibody was obtained by immunization of A4GNT recombinant protein (Accession Number: XM\_017006544). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107668(A4GNT antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human stomach slide using Catalog No:107668(A4GNT Antibody) at dilution of 1:50



# A4GNT antibody



Catalog Number: 107668

---

Immunohistochemistry of paraffin-embedded

human stomach slide using Catalog  
No:107668(A4GNT Antibody) at dilution of 1:50