

# BBS10 antibody

Catalog Number: 108368

## Product name

BBS10 antibody

## Specificity

Human,Zebrafish; other species not tested.

## Antibody description

BBS10 Rabbit Polyclonal antibody. Positive IF detected in hTERT-RPE1 cells. Positive WB detected in HepG2 cells, A549 cells, MCF7 cells. Observed molecular weight by Western-blot: 81kd

## Preparation

This antibody was obtained by immunization of BBS10 recombinant protein (Accession Number: BC013795). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



IF result (nuclear and the base of cilia stain) of anti-BBS10 (Catalog No:108368; 1:50) with hTERT-RPE1 cell (PFA fixed) by Dr. Moshe Kim.



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108368(BBS10 antibody) at dilution of 1:800