

PDCD6IP antibody



Catalog Number: 113768

Product name

PDCD6IP antibody

IP: 1:200-1:2000

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

PDCD6IP Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, mouse liver tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 96kd




Preparation

This antibody was obtained by immunization of PDCD6IP recombinant protein (Accession Number: NM_013374). Purification method: Antigen affinity purified.

IP Result of anti-ALIX; AIP1 (IP:Catalog No:113768, 3ug; Detection:Catalog No:113768 1:500) with Jurkat cells lysate 4000ug.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF, IP, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunofluorescent analysis of HeLa cells, using PDCD6IP antibody Catalog No:113768 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

Jurkat cells were subjected to SDS PAGE followed

PDCD6IP antibody



Catalog Number: 113768

by western blot with Catalog No:113768(ALIX; AIP1 Antibody) at dilution of 1:1000